

Advance 320

CURING LIGHT



- Preset Curing times from 10 to 60 seconds or continuous cure
- LED timed use display in seconds after each cure
- Consistent output and curing strength
- Comes with 8mm light guide and protective eye shield
- Quite fan running-long working time without overheating
- Powerful 75 watt bulb

LIGHTING GUIDES OPTION




101

13mm straight for large facial surface



102

13mm x 80° for posterior occlusal curing



103

8mm x 60° for general curing



104

3mm x 90° for layer curing of posterior restoration